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GB, IRL

Attention: The manual of the offered device with ATEX approval is available in German, English and French (DVR 3: also Dutch). The user is only permitted to operate the device if he understands one of these languages or a technically correct translation of the complete manual. The manual must be completely read and understood before operation of the device and all required measures must be applied. The device has an ATEX approval according to EC-directives and standards. Users outside of the European Community must check national regulations for hazardous areas.

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INSTRUCTIONS MANUAL: EXPLOSION PROOF INDUCTION MOTORS

Contents

INTRODUCTION	2
GENERAL CHARACTERISTICS	2
MOTOR MARKINGS.....	3
USE OF MOTORS	4
INCOMING INSPECTION	4
STORAGE	5
SAFETY CONSIDERATIONS	5
MOUNTING AND START UP	5
MOTOR PROTECTION.....	7
OPERATING CONDITIONS	7
BALL BEARINGS	7
MAINTENANCE.....	7
REPAIR.....	7
SPARE PARTS.....	8
CONDITIONS OF USE.....	9
APPLICABLE STANDARDS.....	10
MANUFACTURER INFORMATION:	10

ELNOR MOTORS NV

INTRODUCTION

An explosion-proof, and in particular a flameproof motor, is one that is inherently safe when working in potentially explosive atmospheres.

This means that its construction and assembly are designed so that ingress of explosive gases or dusts is rendered harmless for the environment.

In the event of accidental ignition, provoking an explosion inside the motor housing, the motor shell and brackets are strong enough to resist the internal pressure created by this explosion.

Moreover, assembly joints and shaft clearance must be long enough and offer minimum play, so that hot gases or dusts resulting from an explosion will be cooled sufficiently to create no possibility of an explosion outside the motor.

The outside temperature of the motor enclosure must also be kept well below the flash point of the gases or dusts present in the hazardous atmosphere, both in normal and abnormal working conditions.

Finally, it is the manufacturer's concern that components inside the motor should be able to withstand an occasional explosion without being adversely affected.

⚠ As explosion proof motors are constructed with special materials and technologies that comply with the legal regulations concerning potentially explosive atmospheres, a wrong connection or a minor modification of the motor destroys the compliance with security regulations. The rules concerning explosion proof apparatus must be observed unconditionally. Note that our motors are approved for a specific group of hazardous areas and temperature classes.

If you are unclear about information's given in this documentation, or if you require further information, please contact Elnor Motors.

GENERAL CHARACTERISTICS

Flame proof induction motor for explosive atmosphere.

Motor prefix	Motor Type	Max. Voltage (VAC)	Frequency (Hz)	Max. power (kW)
BA	3xy***	690	50,60 or 50/60	2,2
BAX	2xy***	400	50,60 or 50/60	0,35
or	1xy***	415	50,60 or 50/60	1
BAV	1xy***H2	415	50,60 or 50/60	1
	73xy***	690	50,60 or 50/60	0,55

BA – naturally cooled explosion proof motor;

BAV – fan cooled explosion proof motor;

BAX - explosion proof motor cooled by application;

xy – length of the stator stack, if length less than 100 mm;

*** - TR – three phase motor;

- C – single phase motor with start capacitor;

- E – single phase motor with permanent run-capacitor;

- CP – single phase motors with start and permanent-run capacitor;

Outside dimensions:

Since our motors are custom made, the outside dimensions, the motor foot/plate, the shaft and the motor face are particular to each customer, please refer to the general drawing attached to your offer.

ELNOR MOTORS NV

MOTOR MARKINGS

Product Marking

The equipment will be marked with the information shown below. Actual layout may vary depending on whether the equipment is single marked or multi marked (ie. ATEX, UKEX, IECEx, INMETRO).

ELNOR MOTORS N.V. B-3150 Haacht Belgium	Company name and address
AB xx ATEX yyy X	AB – name of the notified body issuing the EU type certificate; xx – year of issue of certificate; yyyy- number of EU type certificate; X – special conditions of use;
IECEx AB x x.yyyy X	AB – name of the certified body issuing the IECEx certificate; xx – year of issue of certificate; yyyy- number of Ex certificate; X – special conditions of use;
AB xx UKEX yyyy X	AB – name of the approved body issuing the UKEX certificate; xx – year of issue of certificate; yyyy- number of the UKEX certificate; X – special conditions of use;
II2 G	Equipment group II – non-mining; Equipment category 2 – high protection (see table1); Ignition source risk G –gas;
Exdb IIB	Protection concept Ex d – flameproof enclosure Gas group IIB – Ethylene (for a complete list of gases in this category please refer to standard IEC60079-20-1)
Exdb IIC	Protection concept Ex d – flameproof enclosure Gas group IIC – typical gases are hydrogen and acetylene (for a complete list of gases in this category please refer to standard IEC60079-20-1)
IIB+H2	Gas group IIB+H2 (for a complete list of gases in this category please refer to standard IEC60079-20-1)
Tx	Maximum surface temperature T3, T4, T5 or T6 (see table2)
Gb	Equipment protection level Gb – zone 1(see table1)
Tamb –xx°C to +xx°C	Ambient temperature from minimum -40°C to maximum +70°C, depending on the motor type and protection type.
xxV	Nominal voltage
RPM xxxx	Nominal speed
xx A	Nominal current
P = xxx W	Nominal power
Sx	Service class
CL F	Thermal class
N° JJWW/NFNF xxx	Serial number where: JJ – last 2 digits of the year; WW – number of the week; NFNF – fabrication number; xxx – batch order number;
Customer Ref.	Customer reference if needed;

USE OF MOTORS

⚠ The ELNOR motor is intended for industrial use only. The motor should not be used until the machine in which it is to be installed had met all the necessary approval and regulations. Before installation, the user must check if the group and protection class, marked on the motor label, corresponds with the requested conditions.

ZONES AND EQUIPMENT CATEGORIES				
Zones ¹⁾ ²⁾	Broad definitions of the zones (for guidance only) IEC/EN 60079-10-1	Equipment category 2014/34/EU	Equipment protection level IEC/EN 60079-0	Equipment integrity requirements
1	An area in which it is expected that an explosive gas atmosphere will occur occasionally during normal operation.	2	Gb	Equipment must be safe under normal operation, expected malfunction.

Table 1 zones and equipment categories

⚠ Under no circumstances may Elnor motors be used in zone0

1) - Motors of- Zone 1 can also be used in Zone 2

2) - Motors which are certified for gas must not be used in dust or in hybrid mixtures! Hybrid mixtures: when explosive gas and dust atmospheres occur simultaneously.

TEMPERATURE CLASS	
T CLASS	MAXIMUM SURFACE TEMPERATURE
T1	450°C
T2	300°C
T3	200°C
T4	135°C
T5	100°C
T6	85°C

Table 2 - Temperature class

Temperature classification of the motor depends on the internal thermal protection and it is given on the name plate.

⚠ Under no circumstances should a lower temperature class be used to replace a higher temperature class motor (ie. T3 should not replace a T6).

INCOMING INSPECTION

When receiving the motors, check if any damage has occurred to the packaging of the motors during transportation. If so, immediately notify the carrier, the insurance company and supplier. Failure to do so may invalidate any possible claims.

Inspect the following after packaging has been removed:

- Check for any damage to the motor;
- Check nameplate data to ensure correct motor has been received;
- Check if dimensions and colour are in compliance with your specifications;
- Rotate the shaft by hand to make sure it is turning freely;
- Make sure the motor was not exposed to excessive dirt and moisture during transportation and storage.

In the unlikely event of a claim, please contact our Sales Office.

⚠ Do not use the motor when damaged.

STORAGE

The motor should be stored in a clean, dry and vibrations free environment. If the motor is to be stored for a long period of time, the shaft protected with anti-rust coating should be checked and "touched up" if necessary.

SAFETY CONSIDERATIONS

 Electric motors are source of electrical hazards and consequently their improper use may cause injury and/or damage to people, animals and objects.

Only properly qualified personnel are permitted to perform such tasks as transport, assembly, commissioning and maintenance. The qualified personnel must know and observe all relevant safety standards.

The motor may not be opened in the presence of potentially explosive atmospheres, or when energized.

MOUNTING AND START UP

Electric data of the motor, valid for the rated power, can be found on the nameplate.


It is the responsibility of the user to ensure that the motor nameplate voltage and frequency are the same as the mains supply.

 The user has to ensure that the earth connection is properly made. The grounding connectors are provided inside the terminal box and on the motor frame. According to IEC 60079-0, the grounding cable must have a cross section area of at least 4 mm². Always ensure compliance with local regulations and demands.

  Connection to the mains supply should only be executed by a recognised professional, respecting local demands and regulations.

If motors are supplied without cable the user must ensure that the diameter of the cable is adequate for the voltage and current involved. (See nameplate data of motor)

The connection of the wires must be done as described in the connection diagram that is delivered with the motor. The connection diagram is usually located inside the cover of the terminal box.

The connection diagram is structured as follows: (see example)

The leads can be coloured, numbered or both. Each line represents one lead. The given connection diagram is an example only.

An indication of the tightening torques for the electric connections is give in table bellow:

Thread size	M4	M5	M6
Torque (Nm)	1,5	2	3

Table 3. Tightening torques for electrical connections

The motor must be protected against overloads and short-circuits, through the use of a fuse that corresponds to the current draw as indicated on the nameplate.

ELNOR MOTORS NV

Since our motors are custom made, speed and torque curve are also specific for each motor type and they should be requested from the manufacturer.

The manufacturer cannot be held responsible for damage due to incorrect installation or use.

Free circulation of fresh air for cooling the motor must be guaranteed

The maximum ambient temperature must not exceed the given value that is indicated on the nameplate and in the motor certificate.

To guarantee the explosion proof protection, certified explosion proof cable glands with sealing joints adapted to the cable diameters, must be used in compliance with IEC 60079-0 and IEC 60079-1.

If cable glands are provided with the motor, the cable diameter range as indicated in the attached document, has to be respected. The supplied cable glands are suitable for non-shielded cables only. The cable must be chosen to withstand the maximum cable temperature, as indicated on the nameplate.

⚠ THE MOTOR WILL LOSE ITS EXPLOSION PROOF CHARACTERISTICS IF THE CABLE AND GLANDS ARE NOT FITTED CORRECTLY!

Note: If a spare opening for a cable gland isn't used, only the use of cable glands or blanking elements rated Ex db IIB Tx Gb - where Tx is the temperature class of the motor, as given on the name plate - suitable for use in an ambient of $-40^{\circ}\text{C} \leq T_{\text{amb}} \leq +70^{\circ}\text{C}$, minimum, in compliance with IEC 60079-0 and IEC 60079-1, retains the explosion proof characteristics.

Particular attention is to be given to the fixing torque of bolts and screws. The motor fixing bolts are factory mounted with the correct torque and should not be touched.

It may be necessary to open the terminal box for connecting the motor. When re-mounting the cover, only use the original screws. The fastening torque for the bolts is given in the table below:

Bolt thread size	M5	M6	M8
Torque (Nm)	4 - 5	7 - 9	12 - 14

Table 4. Fastening torque for bolts

- ⚠ Firstly tighten all bolts by hand before applying the full torque. If 4 bolts are present, firstly tighten two diametrically opposed bolts and then the other two.
- ⚠ Motors without a terminal box as all motors in range BA(X)2xx and some motors in range BA(X, V)3xx may not be opened. If opened motors will lose their explosion proof characteristics.
- ⚠ When the threaded blind holes provided on the motor flange are used, only appropriate bolts or screws should be used. The size of the thread engagement, diameter and maximum length, is provided on the customer specific drawing. Under no circumstances should a longer thread engagement be used.

Using a pulley or a coupling device:

Apply only well-balanced discs and couplings. Check the alignment after assembling
Take care while mounting of the pulley on the shaft: Hammering on the shaft will damage the motor bearings. Pulleys should only be mounted with appropriate mounting tools that do not stress the motor bearings.

Following these recommendations will greatly increase the expected motor bearing life.

Safety regulations in respect of guarding couplings and belts from being accidentally touched must be applied.

ELNOR MOTORS NV

Motor type	Diameter (mm)	Length for F_R (mm)	$F_R(N)$ - max		$F_A (N)$ -max
			Cast iron end shield	Aluminium end shield	
BA73xx	11	30		NA	100
BA1xx	14	30	600	NA	150
BA2xx	14	11	NA	1000	140
BA3xx	20	100	600	300	200

Table 5. Permissible forces on shaft

- ⚠ F_R is computed assuming no axial force is present and in normal conditions of use;
- ⚠ F_A is computed assuming no radial force is present;
- ⚠ Since our motors are custom made it is impossible to list in the manual all de different shaft/end shield configurations. If your configuration is not present in the above table please ask the motor manufacturer.

MOTOR PROTECTION

Each motor is provided with an internal thermal protection. If the motor is supplied with a thermal protection, Thermik S06.135 or similar, that is brought out separately, the user must connect the protection in a way that the motor is switched off when the protector switches and doesn't start until an operator starts it manually.

OPERATING CONDITIONS

In general our electric motors are designed for operation up to an altitude of 1000m above sea level and an ambient temperature between -20°C and +40°C. Any variation is stated on the nameplate. Voltage variation may not exceed $\pm 10\%$. Motors should be kept dust free to allow heat exchange to occur.

BALL BEARINGS

The life of the ball bearings and radial and axial forces are motor specific, contact the motor manufacturer for information over life expectancy of ball bearing for your particular application. The maximum operating temperature of the ball bearings is 150C.

MAINTENANCE

Re-greasing of the motor bearings is not necessary, the ball bearing are lubricated for life. Dust build-up on the motor must be avoided.

If the noise level is enhanced, switch off the motor immediately and inspect. Replacement of the bearings must only be performed at the manufacturer.

- ⚠ Unless agreed in writing maintenance action are to be done only by the manufacturer.

REPAIR

The explosion proof motor may not be opened while the motor is energised.

- ⚠ Unless agreed in writing repairs are to be done only by the manufacturer.

ELNOR MOTORS NV

Warning:

 The flameproof joints are not intended to be repaired or reworked by the user.

 The manufacturer is not responsible for motors that are modified or opened without his written agreement.

END OF LIFE

Dismantling of the motor should be done only by qualified personnel following the local safety regulations in place.

Motor contains iron, aluminium, copper and plastic parts which should be recycled by an approved recycling company at the end of life.

SPARE PARTS

For each order of spare parts, the type and number of the motor must be given; these data are indicated on the nameplate.

CONDITIONS OF USE

Motor prefix Motor Type Specific conditions

3xy***

ATEX and UKEX

1. An appropriate connection of the free end of the cable shall be foreseen in case the motor is provided with a permanently connected cable.
2. If operated with a frequency drive, or when a temperature class better than a T3 has to be respected a thermal protection device has to be installed. The following table details the maximum set points of the thermal protection device per temperature class:

Temperature class	Set point windings	Set point end shields/bearings
T3	150	135
T4	135	120

IECEX

1. The flameproof joints are not intended to be repaired.
2. The property class of the special fasteners used to fasten the end shields onto the stator housing

Property class	Part of the special fastener
4.6	Hexagon cap screws
5	Hexagon nut (if used)
8.8	Socket cap screws

3. Model BAV : The motor may be equipped with one or two end shaft(s) or with an external fan. For the fan cooled motors in BAV construction, the motor fan is intended for cooling the motor.
4. Model BAX : When the motor is used to drive a fan, the new assembly, fan and motor unit, requires separate certification.

BA
BAX
BAV

2xy***

ATEX UKEX IECEX

1. Temperature classification of motors is as follows:

Thermal protection location	Switch off thermal protection in C	Temperature at the entry point cable in C	Temperature Class
Winding	135	75	T4
	150	84	T3

2. All safety devices shall operate independently of any measurement or control devices required for operation.
3. No modification must be made to the flamepaths of the enclosure without prior consultation with the manufacturer.

1xy***

1. the property class of the special fasteners used to fasten the end shields onto the stator housing must be at least:8,8;
2. BAX - When the motor is used to drive a fan, the new assembly, fan and motor unit, requires separate certification;
3. Optional push buttons can be installed in the wall of the enclosure;

1xy***H2

1. the property class of the special fasteners used to fasten the end shields onto the stator housing must be at least 12.9;
2. BAX - When the motor is used to drive a fan, the new assembly, fan and motor unit, requires separate certification;

73xy***

1. An appropriate connection of the free end of the cable shall be foreseen when the motor is provided with a cable permanently connected to it;
2. The flamepaths are specified in the manufacturer drawing;

ELNOR MOTORS NV

APPLICABLE STANDARDS

IECEX/IEC	ATEX	UKCA	
IEC 60204-1:2016	EN 60204-1:2018	BS EN 60204-1:2018	Safety of machinery – Electrical equipment of machines – Part 1: General requirements
IEC 60034-1:2017	EN 60034-1:2010	BS EN 60034-1:2010	Rotating electrical machines – Part 1: Rating and performance
IEC 60079-0:2017	EN IEC 60079-0:2018	BS EN IEC 60079-0:2018	Explosive Atmospheres – Part 0: Equipment – General requirements
IEC 60079-1:2014	EN 60079-1:2014	BS EN 60079-1:2014	Explosive Atmospheres – Part 1: Equipment – Equipment protection by flameproof enclosures "d"
	EN 50581:2012	BS EN 50581:2012	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances
	Directives 2014/34/EU		European ATEX Product Directive
		SI 2016 No. 1107 amended by SI2019 No.696	HEALTH AND SAFETY The Equipment and Protective Systems Intended for Use in Potentially Explosive Atmospheres

MANUFACTURER INFORMATION:

ELNOR MOTORS NV
De Costerstraat 45
B-3150 HAACHT
BELGIUM

Telephone: +32 16 60 13 94
Fax: +32 16 60 64 40
E-mail: info@elnor.be
Internet: <http://www.elnor.be>

This document has been made with the greatest care and is correct at the time of printing. It is meant as a general guideline and it is not meant to be complete. Always check with relevant regulations and guidelines.

APRAGAZ
a.s.b.l.
Your inspecting authority

- (1) **EU TYPE EXAMINATION CERTIFICATE**
- (2) **Equipment or Protective system intended for use
in potentially explosive atmospheres
Directive 2014/34/EU**
- (3) EU type examination certificate number: **APRAGAZ 20ATEX0209 X**
- (4) Equipment or protective system: Three-phase and single-phase asynchronous motors
Single-phase: BA(V or X) 3xy C or E or CP
Three-phase: BA(V or X) 3xy TR
- (5) Applicant - Manufacturer: Elnor Motors NV Belgium
- (6) Address: De Costerstraat 45
3150 Haacht (Wespelaar)
Belgium
- (7) This equipment or protective system and any acceptable variation thereto is specified in the schedule to this certificate and the documents therein referred to.
- (8) APRAGAZ, notified body n° 0029 in accordance with article 17 of the directive 2014/34/EU of 26 February 2014, certifies that this equipment or protective system has been found in compliance with the essential health and safety requirements relating to the design and construction of equipment and protective systems intended for use in potentially explosive atmospheres, given in annexe II of the directive.
- The examination and test results are recorded in the confidential file n° 17566/006.
- (9) The conformity of the essential health and safety requirements has been assured by compliance with:
EN 60079-0: 2018
EN 60079-1: 2014
- (10) If the sign "X" is placed after the certificate number, it indicates that the equipment or protective system is subject to special conditions for safe use specified in the schedule to this certificate.
- (11) This EU TYPE EXAMINATION CERTIFICATE relates only to the design and construction of this specified equipment or protective system in accordance with the directive 2014/34/EU. Further requirements of the directive apply to the manufacturing and supply of this equipment or protective system. These are not covered in this certificate.
- (12) The marking of the equipment or protective system shall include the following:
⊕ II 2 G Ex db IIB T4...T3 Gb

Brussels, 07 December 2020

APRAGAZ
Chaussée de Vilvorde 156
1120 Bruxelles
tel 02/264 03 60 fax 02/268 89 58
mail : info@apragaz.com

Ir. B. Nève
General Manager

This certificate can only be reproduced in his entirety, schedule included, without any modification

Rev. 0

1/3

R.A. 001 rev2 - 08-2020

APRAGAZ
a.s.b.l.
Your inspecting authority

- (13) **SCHEDULE**
- (14) **EU TYPE EXAMINATION CERTIFICATE N° APRAGAZ 20ATEX0209 X**
- (15) Description of the equipment or protective system

Three-phase and single-phase asynchronous motors:
Single-phase: BA(V or X) 3xy C or E or CP
Three-phase: BA(V or X) 3xy TR

BA: a naturally cooled explosion proof motor
BAV: a fan cooled explosion proof motor
BAX: a motor cooled by the application

The motor consists of a flameproof enclosure which is made of an aluminium or a steel sheet metal carcass and has aluminium or cast iron end shields.

The fastenings are made by hex cap screw of property class 4.6, nuts of property class 5.0 and socket cap screws of property class 8.8.

The motor may be equipped with two ends shafts and can be electrically directly fed by a single cable or through a connection in the terminal box.

Temperature classes

The motors on S1 duty service cooled by the process have the temperature class T3. The motors have a temperature class T4 if they are equipped with a thermal protection suitable to the temperature class. The motors with natural cooling, the motors energized by converters and the motors on S2 to S9 duty service are provided with a direct thermal protection or a PTC resistance suitable to the temperature class.

Electrical characteristics

Power range: from 0,25 kW to 2,5 kW
Voltage supply: up to 690 Volt – For VAC motors: 50/60 Hz
Insulation class: F

Routine tests

The manufacturer shall make the routine verifications and tests necessary to ensure that the electrical apparatus produced complies with the specifications submitted to the testing station together with the prototype or sample.

Every part of the motor legible for routine testing whether it's made from aluminium, iron cast or steel that is tested at 4*Pref is exempted from overpressure pressure testing, otherwise the routine tests at 11,01 bar need to be executed.

The control of the weld of the steel carcass is always to be tested at 11,01 bar.

A motor of every fabrication number is additionally tested on the test bench. This to verify the maximum surface temperatures on different parts on the motor for the specific type of motor.

Cables resisting up to 120°C should be selected for a T3 and cables resisting up to 100°C should be selected for a T4 motor.

This certificate can only be reproduced in his entirety, schedule included, without any modification

Rev. 0

2/3

R.A. 001 rev2 – 08-2020

APRAGAZ
a.s.b.l.
Your inspecting authority

- (13) **SCHEDULE**
- (14) **EU TYPE EXAMINATION CERTIFICATE N° APRAGAZ 20ATEX0209 X**

(15) Marking:

Elnor Motors NV Belgium
De Costerstraat 45
3150 Haacht (Wespelaar)
Belgium

Year CE + Identification n° of NOBO production control phase
Type: BA(X/V) 3xy
Serial number: 650xxxxx
Manufacturer date, fabrication number and batch number : JJWW/NFNF xxx

II 2 G Ex db IIB T4...T3 Gb
APRAGAZ 20ATEX0209 X

Warning:

Do not open under voltage.

Eventual prescriptions

Choice of cable needs to be compatible with a temperature of X °C.
Ambient temperature range: -40°C up to +70°C.
An adequate warning is to be printed when the motor is delivered without cable glands, plugs or cable.

Note:

This X °C is in function of the temperature class of the motor. X=120°C for a T3 motor and X=100°C for a T4 motor.

Depending on the foreseen duty class of the motor the temperature range can be narrowed.

- (16) The examination and test results are recorded in confidential report n° U194718.

(17) Special conditions for safe use:

An appropriate connection of the free end of the cable shall be foreseen in case the motor is provided with a permanently connected cable.

If operated with a frequency drive, or when a temperature class better than a T3 has to be respected a thermal protection device has to be installed directly in the motor. The following table details the maximum set points of the direct thermal protection device per temperature class:

Temperature class	Set point windings in °C	Set point end shields/bearings in C°
T3	150	135
T4	135	120

(18) Essential health and safety requirements:

The essential health and safety requirements are covered by the conformity to European standards listed at (9).

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Hereby the manufacturer declares that the device is in conformity with the directives:

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The technical documentation in accordance with directive 2014/34/EU (ATEX) is deposited on behalf of VACUUBRAND GMBH + CO KG with the relevant authorized depository:

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Hereby the manufacturer declares that the incomplete machinery is in conformity with the directive:

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The technical documentation in accordance with annex VII B has been compiled. The manufacturer undertakes to submit the technical documentation relating to the incomplete machine to the relevant national authorities as paper mold on request.

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This incomplete machine must not be put into service until the machinery into which it is to be incorporated has been declared in conformity with the provisions of the EC Machinery Directive, the harmonized standards, European standards, or the relevant national standards.

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