VENTILE – VALVES – VANNES

VB M
VB M-B

Hinweise zu Installation und Betrieb
Notes concerning installation and use
Notes concernant l'installation et l'utilisation

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NOTICE

Safety information!

Remove all packing material from the packing box. Remove the product from its packing-box and retain all packaging for inspection until the equipment is inspected and tested. Inspect the equipment promptly and carefully. If the equipment is damaged, notify the supplier and the carrier in writing within three days. State the item number of the product together with the order number and the supplier's invoice number.

Use the equipment **for the intended use only**.

Ensure correct direction of gas flow (see arrow on valve housing).

Comply with all applicable and relevant safety requirements (regulations and guidelines). **Implement the required actions and adopt suitable safety measures.**

Ensure that **maintenance** is done only by suitably trained and supervised technicians. As laid down in the statutory regulations (occupational, health and safety regulations and regulations for environmental protection), equipment returned to the manufacturer can be repaired only under certain conditions.

CAUTION

- Check that mains voltage and current conform with the equipment (see rating plate).
- Check that the system is mechanically stable and that all fittings are secure.
- Avoid overpressure in case inert gas is connected. Provide pressure limitation if necessary.
- Ensure that the materials of the wetted parts are compatible with the processed substances, see "Technical data".



WARNING

- ⚠ **Do not use the equipment if it is damaged.**
- ⚠ Do not permit any **uncontrolled pressurizing. Risk of bursting!**
- ⚠ Pay attention to the maximum permissible ambient temperature.
- ⚠ **The electric wiring must be carried out only by suitably trained and supervised personnel!**
- ⚠ Use only **OEM spare parts and accessories**. Otherwise safety and performance of the equipment might be reduced. Ensure that maintenance is done only by suitably trained and supervised technicians.
- ⚠ Make sure that air/gas inlet through the valve **never leads to hazardous, explosive or otherwise dangerous mixtures**. If in doubt, use inert gas.

DANGER

- ➔ Air and processed media might react in the system and form hazardous or explosive mixtures, when you use air rather than inert gas at the gas inlet. This constitutes a risk of significant damage to equipment and/or facilities, a risk of personal injury or even loss of life.

The valves **VB M KF 16 (20666817)** and **VB M G 1/4 (20666818)** are shipped with a cable with a three pole male connector.

The valve **VB M KF 10 (20666816)** has to be suitably wired by the customer:

- ➔ Loosen the screw and remove the cable head from the solenoid operated valve block.
- ➔ Take out the inner part of the cable head.
- ➔ Screw a bifilar cable to the terminals at the inner part of the cable head: terminal 1 (ground) and terminal 2 (24V=).
- ➔ Run the connecting cable through the cable head and through the cable fitting.
- ➔ Screw the cable head with its inner part and with the seal to the solenoid operated valve block.
- ➔ Wire the connection cable with a suitable plug and / or control device.

Technical data

Type	20666816 VB M KF 10	20674217 VB M-B KF 16	20666817 VB M KF 16	20666818 VB M G1/4
Maximum permissible range of electrical supply ($\pm 10\%$)	24 V DC			
Power draw (operation)	W	4		
Electrical connection	terminal inside cable head	connector for VACUU•BUS system	male connector 3-pole, cable length 2.5 m	
Rated duty	continuous operation			
Status of valve when currentless	closed			
Duty ratio	100%			
Switching time		5-18 (pressure-dependent)		
Opening	ms	approx. 8		
Closing	ms			
Nominal diameter	mm	2.4		
Gas throughput	l/min	140		
Maximum permissible pressure	bar	1.5 (absolute)		
Connection (towards application)	small flange DN 10	small flange DN 16	thread G1/4	
Connection (gas inlet)	connection for hose 6/8 mm			
Permissible ambient temperature	°C	0 to +40		
Maximum permissible temperature of gaseous media	°C	-10 to +40		
Dimensions L x W x H approx.	mm	80 x 62 x 44		
Length of connecting cable	m	-	2.5	
Weight approx.	kg	0.25	0.30	

Attention: The internal venting valve of the vacuum controller CVC 3000 is without function if an external venting valve VB M-B is connected!

Components	Wetted materials
Seals	FPM
Housing, seat, inner parts of valve	stainless steel

We reserve the right for technical modifications without prior notice!

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VACUUBRAND®

DECLARATION OF CONFORMITY – China RoHS 2

VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB/T 26572 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(+VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmentally Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmentally Friendly Use Period for VACUUBRAND products is 40 years.

MATERIAL CONTENT DECLARATION FOR VACUUBRAND PRODUCTS						
有毒有害物质或元素 Hazardous substances						
部件名称 Part name	铅 Pb	汞 Hg	镉 Cd	六价铬 Cr(+VI)	多溴联苯 PBB	多溴二苯醚 PBDE
包装 Packaging	0	0	0	0	0	0
塑料外壳 / 组件 Plastic housing / parts	0	0	0	0	0	0
真空油 Vacuum oil	0	0	0	0	0	0
电池 Battery	0	0	0	0	0	0
玻璃 Glass	X	0	0	0	0	0
电子电气组件 Electrical and electronic parts	X	0	0	0	0	0
控制器 / 测量设备 Controller / measuring device	X	0	0	0	0	0
金属外壳 / 组件 Metal housing / parts	X	0	0	0	0	0
电机 Motor	X	0	0	0	0	0
配件 Accessories	X	0	0	0	0	0
此表格是按照SJ/T 11364-2014中规定所制定的。 This table is created according to SJ/T 11364-2014.						

VACUUBRAND®

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Batteries, glassware and accessories might not be content of the enclosed device and may have its own EFUP-marking and/or might be maintaining parts with changing EFUP-marking.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).

With these products and appliances in particular, please note the EFUP labeled on these products.

VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 06 September 2022

(Dr. Constantin Schöler)
Managing Director



ppa. (Jens Kaibel)
Technical Director

VACUUBRAND GMBH + CO KG

Alfred-Zippe-Str. 4
97877 Wertheim
Germany

Tel.: +49 9342 808-0
Fax: +49 9342 808-5555
E-Mail: info@vacuubrand.com
Web: www.vacuubrand.com

[VACUUBRAND > Support > Manuals](#)

Manufacturer:

VACUUBRAND GMBH + CO KG
Alfred-Zippe-Str. 4
97877 Wertheim
GERMANY

Phone:

Head office +49 9342 808-0
Sales +49 9342 808-5550
Service +49 9342 808-5660

Fax: +49 9342 808-5555

Email: info@vacuubrand.com

Web: www.vacuubrand.com