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
**MD 4C EX VARIO**  
**MV 10C EX VARIO**  
**MD 4C EX VARIO + AK + EK**  
**MV 10C EX VARIO + AK + EK**

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2. Serial no.: .....

3. Reason for return / malfunction: .....

4. Has the device been used in a copper process step (e.g., semiconductor production):  yes  no

5. Substances (gases, liquids, solids, biological material, e. g. bacteria, viruses) in contact with the device / which have been pumped:

.....  
 .....  
 .....

6. Risk level of the used biological material:  none  1  2\*  3\*\*  4\*\*

\* Contact the VACUUBRAND service absolutely before dispatching the device.

\*\* Devices which have been in contact with biological substances of risk level 3 or 4 cannot be checked, maintained or repaired. Also decontaminated devices must not returned to VACUUBRAND due to a residual risk.

7. Radioactive contamination:  yes  no

8. Prior to return to the factory the device has been decontaminated:  yes  no

Description of the decontamination method and the test / verification procedure:

.....  
 .....

9. All parts of the device are free of hazardous, harmful substances:  yes  no

10. Protective measures required for service staff:

.....

11. If the paint is damaged, we wish a repaint or a replacement of parts for reason of appearance (repaint and replacement at customer's expense):  yes  no

### 12. Legally binding declaration

We assure for the returned device that all substances, which have been in contact with the device are listed in section 5 and that the information is complete and that we have not withheld any information. We declare that all measures - where applicable - have been taken listed in section "Repair - Maintenance - Return - Calibration". By our signature below, we acknowledge that we accept liability for any damage caused by providing incomplete or incorrect information and that we shall indemnify VACUUBRAND from any claims as regards damages from third parties. We are aware that as expressed in § 823 BGB (Public Law Code of Germany) we are directly liable for injuries or damages suffered by third parties, particularly VACUUBRAND employees occupied with handling/repairing the product. Shipping of the device must take place according to regulations.

Name: ..... Signature: .....

Job title: ..... Company's seal:

Date: .....

Release for repair grant by VACUUBRAND (date / signature):.....  
 Protective measures:  Protective gloves, safety goggles  Hood  External cleaning

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e-mail :  
direction@issep.be  
site web :  
http://www.issep.be

(1) **EC TYPE EXAMINATION CERTIFICATE**

(2) **Equipment or protective system intended for use  
in potentially explosive atmospheres  
Directive 94/9/EC**

(3) EC type examination certificate number: **ISSEP09ATEX014X**

(4) Equipment : Asynchronous flameproof motors  
Series BA(V/X) 3xy \* \* \*

(5) Applicant – Manufacturer:  
ELNOR MOTORS N.V.

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(7) This equipment or protective system and any acceptable variation thereto is specified in the schedule to this certificate and the documents therein referred to.

(8) ISSEP, notified body n° 492 in accordance with article 9 of the Council Directive 94/9/EC of 23 March 1994, certifies that this equipment or protective system has been found to comply with the Essential Health and Safety Requirements relating to the design and construction of equipment and protective systems intended for use in potentially explosive atmospheres given in annex II to the Directive.

The examination and test results are recorded in confidential report no 07087

(9) Compliance with the Essential Health and Safety Requirements has been assured by compliance with:  
EN 60079-0:2006 (IEC 60079-0:2004) IEC 60079-0:2007  
EN 60079-1:2007 (IEC 60079-1: 2007)

(10) If the symbol "X" is placed after the certificate number, it indicates that the equipment or protective system is subject to special conditions for safe use specified in the schedule to this certificate

(11) This EC TYPE EXAMINATION CERTIFICATE relates only to the design, examination and tests of the specified equipment or protective system in accordance to the Directive 94/9/EC. Further requirements of this Directive may apply to the manufacturing process and supply of this equipment or protective system. These are not covered by this certificate.

(12) The marking of the equipment or protective system shall include the following indications:  
 II 2 G Ex d IIB T6 to T3 (according to standard EN 60079-0 : 2006)  
 Ex d IIB T6 to T3 Gb (according to standard IEC 60079-0 : 2007)  
 Colfontaine, the 04.03.2009

INSTITUT SCIENTIFIQUE DE SERVICE PUBLIC  
Zoning A. Schweitzer - B7340 Colfontaine (Wasmes)  
Tél: ++ 32 65 610811 – Fax: ++ 32 65 610808

  
Marcel Lambert,  
Manager.

This certificate may only be reproduced in its entirety and without any change, schedule included

(13)	<b>SCHEDULE</b>
(14)	<b>EC TYPE EXAMINATION CERTIFICATE No ISSeP09ATEX014X</b>
(15)	<p>Description of the equipment or protective system:          Three-phase asynchronous motor :BA(V/X) 3xyTR          Single-phase asynchronous motor :BA(V/X) 3xy C or E or CP          The flameproof enclosure of this motor consists of a carcass made of steel and end shields made of cast iron or aluminium alloy.          The fastenings are made by M8 bolts of property class 4.6 ,M8 internal socket head screws of property class 8.8 and M8 nuts of property class 5          The motor may be equipped with two ends shaft.</p> <p><u>Temperature classes</u>          The motors on S1 duty service cooled by the process have the temperature class T3. These motors have the temperature classes T6 or T5 or T4 if they are equipped with a thermal protection suitable to the temperature class.          The motors with natural cooling &amp; the motors energised by converters and the motors on S2 to S9 duty service are provided with a thermal protection or a PTC thermistance suitable to the temperature class - to see documents from the manufacturer</p> <p><u>Electrical characteristics</u>          Power range : up to 3 kW          Voltages supply : up to 695 VAC - 50 / 60 Hz          Insulation class : F</p> <p><u>Eventual prescriptions</u>          Ambient temperature range: -40°C up to +70°C</p> <p><u>Routine test</u>          Each apparatus shall be submitted to the routine verifications and tests necessary to ensure that it complies with the specification submitted to the testing station with the prototype.          The motors equipped with cast iron end shields are exempted from the overpressure routine test. However a control of the welding of the carcass shall be made.          The motors equipped with aluminium alloy end shields shall be submitted to an overpressure routine test under the value of 1101 kPa (11.01 bar).</p> <p>(16) Report no 07087 dated 2.03.2009 composed in total of 27 pages.          The letter from the manufacturer Reference ES:8.030 of 27 October 2008          The manual of instructions of the manufacturer of 8.11.2007 signed on 10.01.2008 (5 pages)          The data sheets related to the parts of the flameproof enclosure dated 4 January 2008 - Rev of 27 February 2009 signed on 27.02.2009 (4 pages)          The data related to the self-adhesive marking plate - document of 10.05.2007 signed on 4.07.2007 (1 page)          The document related to the minimum radial clearance (k) and the maximum radial clearance (m)of the shaft glands, dated 10.01.2008 (1 page)          The document related to the discharge time of the capacitor dated 25.08.2008 (1 page)</p>
This certificate may only be reproduced in its entirety and without any change, schedule included	

## SCHEDULE

## EC TYPE EXAMINATION CERTIFICATE No ISSeP09ATEX014X

- (16) The drawings  
240000423A of 10.01.2008 Rev.A of 15.09.2008  
507960350E Rev.E of 2.07.2007  
507991425E Rev.E of 7.11.2008  
806123063 of 10.01.2008  
816122880B of 26.11.2001 Rev.B of 24.02.2009  
816122881 of 26.11.2001  
816122883 of 12.12.2001  
816122884A of 14.12.2001 Rev.A of 20.02.2009  
835 C of 23.05.1995
- (17) Special conditions for safe use: symbol X
- An appropriate connection of the free end of the cable shall be foreseen when the motor is provided with a supply cable which is permanently connected to it.
- (18) Essential Health and Safety Requirements: covered by the standards listed under point 9 of the present certificate and by the descriptive documents from the manufacturer.

This certificate may only be reproduced in its entirety and without any change, schedule included

## VARIATION

## EC TYPE EXAMINATION CERTIFICATE N° ISSeP09ATEX014X/1

- (14) Equipment:  
Three-phase asynchronous motors: BA(X)3xy TR  
Single-phase asynchronous motors: BA(X)3xy C or E or CP
- (15) Subject to the variation:  
To allow the mounting of the aluminium alloy cover described on the drawing 517994625 for the version of motor equipped with a terminal box.  
The motors equipped with this cover are exempted from the overpressure routine test.  
However a control of the welding of the carcass shall be made
- Marking: Unchanged
- Eventual prescriptions  
The ambient temperature range: -40°C up to +70°C
- (16) Report n° 09131 of 26.01.2010 (16 pages)  
The data sheets related to the construction elements of the flameproof enclosure dated 15.01.2010 (5 pages)  
The e-mail from the manufacturer of 21 January 2010  
The drawing  
517994625 of 2.11.2009
- (17) Special conditions for safe use: Unchanged
- (18) Essential Health and Safety Requirements: Covered by the standards listed at point 9 of the initial certificate

Colfontaine, 29.01.2010

  
Lambert Marcel,  
DirecteurINSTITUT SCIENTIFIQUE DE SERVICE PUBLIC  
Zoning A. Schweitzer - B 7340 Colfontaine (Wasmès)  
Tel: ++ 32 65 610811 - Fax: ++ 32 65 610808

This document may not be used without the original certificate

## VARIATION

### EC TYPE EXAMINATION CERTIFICATE N° ISSeP09ATEX014/2

- (14) Equipment:  
Three-phase asynchronous motors: BA(X)3xy TR  
Single-phase asynchronous motors: BA(X)3xy C or E or CP
- (15) Subject to the variation:  
To allow the mounting of aluminium alloy end shields described on the drawings 517412215D and 517431005B  
To allow the mounting of an aluminium alloy carcass ( stator housing)  
To attest the update of several drawings  
The conformity assessment has been made on basis of the standard IEC 60079-0 edition 2011.  
  
Marking: Unchanged except the indication of the range of temperature classes in the code: T6...T3  
  
Eventual prescriptions  
The ambient temperature range: -40°C up to +70°C
- (16) Report n° 11122 of 11.07.2012 (15 pages)  
The data sheets related to the construction elements of the flameproof enclosure dated 29.06.2012 (4 pages)  
The e-mail from the manufacturer of 11.07.2012  
The drawings  
107991400G of 3.06.2010  
117412200F of 20.09.2011  
117431000C of 6.09.2011  
117994600H of 19.03.2010  
230000281C of 1.06.2005  
507960325D of 21.12.2006  
507960345D of 23.09.2010  
507991425F of 19.02.2010  
517412215D of 7.07.2011  
517431005B of 7.09.2011  
517992201J of 1.12.2010  
517994601H of 20.10.2009  
816122883A of 5.12.2011  
816122884C of 29.11.2011  
816123090 of 18.02.2010
- (17) Special condition for safe use: Unchanged
- (18) Essential Health and Safety Requirements: Covered by the standards listed at point 9 of the initial certificate

Colfontaine, 17.07.2012

INSTITUT SCIENTIFIQUE DE SERVICE PUBLIC  
Zoning A. Schweitzer, rue de la Platerie  
B-7340 Colfontaine (Wasmes)  
Tél: ++ 32 65 610811 – Fax: ++ 32 65 610808  
e-mail : colfontaine@issep.be

  
F. J. M. Lambert  
Directeur

This document may not be used without the original certificate

## VARIATION

## EC TYPE EXAMINATION CERTIFICATE N° ISSeP09ATEX014X/3

(14) Equipment

- Three-phase asynchronous motors: BA (V/X) 3xy TR.
- Single-phase asynchronous motors: BA (V/X) 3xy C or E or CP.

(15) Subject to the variation

- To allow the mounting of two new endshields described on the drawings 117992200K et 117992900.
- To attest the update of drawings.

Marking: Unchanged.

Eventual prescriptions: Unchanged.

(16) Report N° 12116 completed by 8 documents.(17) Special condition for safe use: Unchanged.(18) Essential Health and Safety Requirements: Covered by the standards listed at point 9 of the initial certificate.

Colfontaine, 03.04.2013.

INSTITUT SCIENTIFIQUE DE SERVICE PUBLIC  
Zoning A. Schweitzer, rue de la Platinerie  
B-7340 Colfontaine (Wasmes)  
Tél: ++ 32 65 610811 – Fax: ++ 32 65 610808  
colfontaine@issep.be

  
Marcel LAMBERT,  
Director.

This document may not be used without the original certificate

## DECLARATION OF CONFORMITY – China RoHS 2

VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(+VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmentally Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmentally Friendly Use Period for VACUUBRAND products is 40 years.

此表格是按照SJ/T 11364-2014中规定所制定的。  
This table is created according to SJ/T 11364-2014.

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Declaration of Conformity – China RoHS 2

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).

With these products and appliances in particular, please note the EFUP labeled on these products.

VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

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Copyright 2020