

VacuHandControl VHC^{pro}

[Link to the website](#)



Product description

The VHC^{pro} is an ergonomic light weight suction hand set, fully autoclavable at 121 °C. It is ready for use with 2.5 m flow through and very flexible silicon tube, adapters for graduated, aspiration and pasteur pipettes, adapter for pipette tips, a wall holder for standby-position (keeps the table top free) and clips with adhesive surface to guide the suction tube into the biosafety cabinet neatly. The mechanism of the VHC^{pro} handset is not in contact with contaminated fluids and enables reliable, leak free working.

Performance features

- light weight and ergonomic for convenient handling
- permanent suction mode for continuous flow
- flow through media tube for leak free working
- autoclavable at 121 °C and 2 bars

Ordering data

Product name	Nominal mains voltage mains frequency	Plug type	Articleno.
VHC ^{pro}			20688061

Items supplied

Ergonomic aspiration hand set with tip holder and adapter to accept a variety of common pipettes and tips, level for manual suction control, rotary knob for permanent aspiration, with 2.5 m tubing.

Accessories

Technical data

Permissible ambient temperature range (operation)	10 - 40 °C
Permissible ambient temperature range (storage)	-10 - 60 °C
Max. outlet pressure (abs.)	1.1 bar
Inlet connection	Suction controller
Outlet connection	--
Protection class IEC 60529	IP 20
NRTL certification	No

Dimension sheet



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com