

VacuHandControl VHC^{pro}

[Enlace al sitio web](#)



Descripción del producto

El VHC^{pro} es un ligero y ergonómico mando de succión, completamente autoclavable a 121 °C. Está preparado para usarse con un tubo muy flexible de silicona de 2.5 m con adaptadores graduados, pipetas de aspiración y pasteur, adaptador para puntas de pipeta, un soporte de pared (mantiene la mesa libre) y clips con superficie adherente para guiar el tubo de succión en el interior de la cabina de bioseguridad cuidadosamente. El mecanismo del mando VHC^{pro} no está en contacto con fluidos contaminados y está realmente libre de fugas.

Rendimiento

- ligero y ergonómico para un conveniente manejo
- succión permanente para un flujo continuo
- tubo libre de fugas
- autoclavable a 121 °C y 2 bares

Información de pedido

Descripción del producto	Tensión nominal frecuencia de red	Tipo de enchufe	N° de artículo
VHC ^{PRO}			20688061

Artículos suministrados

Mango de aspiración ergonómico con soporte de puntas y adaptador para aceptar una variedad de pipetas y puntas comunes, palanca operativa para la aspiración manual y mando giratorio para la aspiración permanente, con 2.5 m de manguera.

Accesorios

Datos técnicos

Temperatura ambiente admisible durante (operación)	10 - 40 °C
Temperatura ambiente admisible durante (almacenaje)	-10 - 60 °C
Max. Presión de salida (abs.)	1.1 bar
Conexión, lado de aspiración	Dispositivo de succión
Conexión, lado de la presión	--
Clase de protección según IEC 60529	IP 20
Certificación NRTL	No

Hoja de dimensiones



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com