

真空传感器

Cerabar PMP21 KF DN 16

网站链接



产品说明

性能特点

- ATEX规格：II 2G Ex ia IIC T4 Gb
- 测量范围1000-1 mbar
- 精度：好于+3 mbar

订购信息

产品名称	额定电源电压 工作频率	电源插头类型	订货号
Cerabar PMP21 KF DN 16	VACUU-BUS DC 24 V DC	ISO4400/M16	20635424

货品配置

该真空传感器用于危险区域，完全配备KF DN 16小法兰和插头电缆组件

配件

技术指标

测量范围 (mbar)	1000 - 1 mbar
测量范围 (torr)	750 - 1 torr
测量原理	Edelstahl-Membran (1.4435), gasartunab., Absolutdruck
测量的不确定性	< +/- 3 mbar
真空接头	KF DN 16 小法兰
的允许环境温度范围 (运行)	10 - 40 °C
的允许环境温度范围 (存放)	-10 - 60 °C
最大介质温度 (持续操作)	40 °C
短期内介质温度的最大值	80 °C
外壳材质	金属外壳
防护等级	IP 65
ATEX规格	II 2G Ex ia IIC T4 Gb
NRTL 标准	否
通讯	4 - 20 mA

尺寸规格表



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com