

操作面板 VACUU·SELECT

[网站链接](#)

产品说明

这只是 VACUU·SELECT 控制器的 VACUU·SELECT 操作单元，用于特殊安装和更换情况。供货范围不包括真空传感器或电源装置。如需全套 VACUU·SELECT 控制器，请订购货号 20700000 的产品。其他型号的控制器请参见紧凑型控制器。用于精细真空控制的 VACUU·SELECT 套件见精细真空控制套件部分。

性能特点

- 只有 VACUU·SELECT 控制单元适用于特殊安装和更换情况
- 无真空传感器
- 不带电源装置
- 如需购买配备齐全的 VACUU·SELECT 控制器，请订购货号 20700000 的产品

订购信息

产品名称	额定电源电压 工作频率	电源插头类型	订货号
VACUU·SELECT	VACUU·BUS DC	不包括	20700040

货品配置

配件

- VACUU·SELECT 用户管理 VACUU·SELECT User Management (20700900)
- 橡胶真空管 DN 8 mm (20686001)

技术指标

的允许环境温度范围（运行）	10 - 40 °C
的允许环境温度范围（存放）	-10 - 60 °C
最大介质温度（持续操作）	40 °C
短期内介质温度的最大值	80 °C
外壳材质	具有良好化学耐腐蚀性的坚固的塑料外壳
防护等级	IP 40
防护等级UL50E	Type 1
依据IEC 60529前面板的防护等级	IP 41
依据UL50E前面板的防护等级	Type 2
尺寸（长×宽×高）	152 x 127 x 41 mm
重量	0.59 kg
最大阻断容量/阻断电流, 24V/VACUU-BUS®	4 A
连接电缆长度	ohne Netzkabel m
NRTL 标准	是
通讯	Ethernet (RJ45): Modbus TCP / USB (Type A): Data transfer, software update, RS-232 (additional adapter required)

尺寸规格表



VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com